




Gynecology & Obstetrics


- **Single Drapes**
- **Laparoscopy Packs**
- **Caesarean Packs**
- **Delivery Packs**
- **Lithotomy Packs**

LAPAROSCOPY DRAPE WITH ATTACHED LEGGING




Product Code #	Description	Size(cm)
43101-000	Laparoscopy Drape with Attached Legging <ul style="list-style-type: none"> • 15x20 cm Perineal Fenestration • Absorbent Laminated Reinforcement • Adhesive • Leggings • 30x30 Abdominal Fenestration 	200x160

LAPAROSCOPY T" DRAPE



Product Code #	Description	Size(cm)
43102-000	Laparoscopy T" Drape <ul style="list-style-type: none"> • 15x7 cm Perineal Fenestration with Flap • 15x20 Abdominal Fenestration • Absorbent Laminated Reinforcement • Tube Holder • Adhesive 	240x220


ABDOMINAL DRAPE



Product Code #	Description	Size(cm)
40001-000	Abdominal Drape <ul style="list-style-type: none"> • Adhesive • Absorbent Laminated Reinforcement 	150x100




LEGGINGS




Product Code #	Description	Size(cm)
CC 142005-002	Leggings Plastic Blue • Latex free elastic closure	140x73

CONE SHAPE LEGGINGS



Product Code #	Description	Size(cm)
CC 142005-004	Cone Shape Leggings	104x110


CONE SHAPE FOOT COVER



Product Code #	Description	Size(cm)
CC 142005-005	Cone Shape Foot Cover • Latex Free Elastic Closure	19x140




LEGGINGS BLUE PLASTIC




Product Code #	Description	Size(cm)
CC 142005-003	Leggings Blue Plastic • Latex Free Elastic Closure	55x38

MATERNITY MAT WITH POUCH



Product Code #	Description	Size(cm)
CR 143038-007	Maternity Mat with Pouch • Under Pad • Laminated Sheet	140x65


MATERNITY MAT WITH FLUID COLLECTION POUCH



Product Code #	Description	Size(cm)
CR 143038-004	Maternity Mat with Fluid Collection Pouch • Under Pad • Fluid Collection Pouch	175x65




LEGGINGS




Product Code #	Description	Size(cm)
CC 142005-000	Leggings	160x100

UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH



Product Code #	Description	Size(cm)
CR 143038-005	Under Buttocks Drape with Fluid Collection Pouch <ul style="list-style-type: none"> • Adhesive • Absorbent Laminated Reinforcement • Fluid Collection Pouch 	75x120


UNDER BUTTOCKS DRAPE WITH ABSORBENT PAD



Product Code #	Description	Size(cm)
CR 143038-006	Under Buttocks Drape with Absorbent Pad <ul style="list-style-type: none"> • Absorbent Pad 	75x120




LITHOTOMY T™ SHEET




Product Code #	Description	Size(cm)
30226-0000	Lithotomy T™ Sheet <ul style="list-style-type: none"> • Adhesive • Absorbent Laminated Reinforcement • Tube Holder • 7x15 Cm Perineal Fenestration 	190x150

C-SECTION DRAPE



Product Code #	Description	Size(cm)
30162-0000	C-Section Drape <ul style="list-style-type: none"> • Adhesive • 20x20 Cm Tear Drop Fenestration • Absorbent Laminated Reinforcement • Tube Holder • Arm Board Covers 	300x160


CAESAREAN DRAPE WITH FLUID COLLECTION POUCH



Product Code #	Description	Size(cm)
30161-0000	Caesarean Drape with Fluid Collection Pouch <ul style="list-style-type: none"> • Hygiene Sheet • Fluid Collection Pouch • Incise • 30 Cm Circular Fenestration 	300x160




LARGE SHEET




Product Code #	Description	Size(cm)
41013-0000	Large Sheet <ul style="list-style-type: none"> • Absorbent Laminated Reinforcement • Adhesive • Tube Holder 	160x200

MEDIUM SHEET




Product Code #	Description	Size(cm)
41014-0000	Medium Sheet <ul style="list-style-type: none"> • Absorbent Laminated Reinforcement • Adhesive • Tube Holder 	160x100

SMALL SHEET



Product Code #	Description	Size(cm)
41015-0000	Small Sheet <ul style="list-style-type: none"> • Absorbent Laminated Reinforcement • Adhesive • Tube Holder 	100x80


EXTRA SMALL SHEET



Product Code #	Description	Size(cm)
41016-0000	Extra Small Sheet <ul style="list-style-type: none"> • Absorbent Laminated Reinforcement • Adhesive 	80x80




SURGICAL SHEET MEDIUM




Product Code #	Description	Size(cm)
40015-0000	Surgical Sheet Medium • Adhesive • Absorbent Laminated Reinforcement	160x120

SURGICAL SHEET




Product Code #	Description	Size(cm)
40017-0000	Surgical Sheet • Adhesive • Absorbent Laminated Reinforcement	80x80

SURGICAL SHEET SMALL



Product Code #	Description	Size(cm)
40016-0000	Surgical Sheet Small • Adhesive • Absorbent Laminated Reinforcement	100x100

SMALL FENESTRATED SHEET




Product Code #	Description	Size(cm)
40011-0000	Small Fenestrated Sheet • 8 Cm Circular Fenestration • Absorbent Laminated Reinforcement	100x80




Laparoscopy Pack I : Product Code - #66481-0001

LAPAROSCOPY DRAPE WITH ATTACHED LEGGING




CAMERA COVER




Item Code	Qty.	Contents	Size(cm)
318-080-100S	1 PC	Back Table Cover Reinforced	100x80
CR 143027-024	3 PCS	Surgical Gown	L
CR 143029-000	3 PCS	Hand Towel	30x30
43101-0000	1 PC	Laparoscopy Drape with Attached Leggings <ul style="list-style-type: none"> • 150x20 Cm Perineal Fenestration • 30x30 Cm Abdominal Fenestration • Absorbent Laminated Reinforcement • Adhesive • Leggings 	200x160
CR 143010-002	1 Pc	Camera Cover Plastic	20x200
CP 141066-000	1 Pc	Carry Bag	

Laparoscopy Pack II: Product Code - #66481-0002

LAPAROSCOPY T" DRAPE



LEGGINGS




Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
CR 143027-024	3 PCS	Surgical Gown	L
CR 143029-000	3 PCS	Hand Towel	30x30
CR 143038-006	1 PC	Under Buttocks Drape with Absorbent Pad <ul style="list-style-type: none"> • Absorbent Pad 	75x120
CR 143029-000	1 PC	Hand Towel	30x30
043-060-040A	2 PCS	Utility Drape <ul style="list-style-type: none"> • Adhesive 	60x40
CC 142005-000	2 PCS	Leggings	160x100
43102-0000	1 PC	Laparoscopy T" Sheet <ul style="list-style-type: none"> • 15x7 Cm Perineal Fenestration with Flap • 15x20 Cm Abdominal Fenestration • Absorbent Laminated Reinforcement • Tube Holder • Adhesive 	240x220
CP 141066-000	1 PC	Carry Bag	



Caesarian Pack I : Product Code - #66451-0001

CAESAREAN DRAPE WITH FLUID COLLECTION POUCH




DAITHERMY BAG
40 Cm



Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
CR 143027-024	1 PC	Surgical Gown	L
CR 143029-000	1 PC	Hand Towel	30x30
CR 143027-030	2 PCS	Surgical Gown Reinforced	L
CR 143029-000	2 PCS	Hand Towel	30x30
CR 143041-006	10 PCS	Gauze Swabs with X-ray Detectable	45x45
CR 143123-000	1 Pc	Connecting Tubing 10" Feet	
CR 143120-000	1 PC	Yenkuear Tip	
CR 143120-000	1 PC	Bulb Syringe	
CR 143104-000	2 PCS	Cord Clamp Lock	
CC 142002-000	1 PC	Diathermy Bag	30x35
		• Adhesive	
CR 143039-000	1 PC	Mayo Stand Cover	80x80
30161-0000	1 PC	Caesarian Drape with Fluid Collection Pouch	
		• Hygiene Sheet	
		• Fluid Collection Pouch	
		• Incise	
		• 30 Cm Circular Fenestration	
CR 143033-002	1 PC	ID Band	Infant
CR 141066-000	1 PC	Carry Bag	

Caesarian Pack II : Product Code - #66451-0002

CAESAREAN DRAPE WITH FLUID COLLECTION POUCH




Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
CR 143027-024	3 PCS	Surgical Gown	L
CR 143029-000	3 PCS	Hand Towel	30x30
CR 143041-012	6 PCS	Mopping Pad	30x30
30161-0000	1 PC	Caesarian Drape with Fluid Collection Pouch	
		• Hygiene Sheet	
		• Fluid Collection Pouch	
		• Incise	
		• 30 Cm Circular Fenestration	
CR 143039-000	1 PC	Mayo Stand Cover	80x80
CP 141066-000	1 PC	Carry Bag	




Caesarian Pack III : Product Code - #66451-0003

FENESTRATED ABSORBENT SHEET




ABSORBENT SHEET




Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
CR 143027-024	3 PCS	Surgical Gown	L
CR 143029-000	3 PCS	Hand Towel	30x30
065-150-100S	2 PCS	Absorbent Sheet	150x100
065-100-075A	2 PCS	Absorbent Sheet	100x75
065-150-225S	1 Pc	Absorbent Fenestrated Sheet • Adhesive	150x225
CR 143243-002	1 PC	ISO Seal Drape Tape Strips	5x20
CR 143041-012	6 PCS	Mopping Pad	30x30
CR 143084-000	1 PC	Under Pad	60x60
CR 143039-000	1 PC	Mayo Stand Cover	80x80
043-160-100S	1 PC	Surgical Sheet	160x100
CR 143005-001	1 PC	Apron	L
043-080-080S	1 PC	Surgical Sheet	80x80
CR 143033-002	1 PC	ID Band Infant Pink	
CR 143014-000	1 PC	Cord Clamp Lock	
CP 141066-000	1 PC	Carry Bag	

Caesarian Pack IV : Product Code - #66451-0004

C-SECTION DRAPE




MAYO STAND COVER





Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
CR 143027-030	3 PCS	Surgical Gown Reinforced	L
CR 143029-000	3 PCS	Hand Towel	30x30
CR 143030-004	1 PC	Mayo Stand Cover	60x140
CR 143041-012	6 PCS	Mopping Pad	30x30
CR 143041-006	2 PC	Gauze Swap X-ray Detectable	45x45
30162-0000	1 PC	C-Section Drape *Absorbent Laminated Reinforcement *Adhesive *20x20 Cm Tear Drop Fenestration *Tube Holder *Arm Board Covers	
CR 143121-000	1 PC	Bulb Syringe	
CR 143014-000	1 PC	Cord Clamp Lock	
065-075-075S	1 PC	Absorbent Sheet	75x75
CR 143033-003	1 PC	ID Band Infant Pink	
CR 143005-001	1 PC	Apron	Large
CP 141066-000	1 PC	Carry Bag	



Delivery Pack I : Product Code - #66551-0001





UNDER BUTTOCKS DRAPE WITH FLUID COLLECTION POUCH



Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
043-110-200S	1 PC	Set Up Cover	110x200
CR 143027-024	1 PC	Surgical Gown	L
CR 143029-000	2 PCS	Hand Towels	30x30
065-075-075S	1 PC	Absorbent Sheet	75x75
CR 143038-005	1 PC	Under Buttocks Drape with Fluid Collection Pouch	120x75
CC1 42005-000	2 PCS	Leggings with 15 Cm Cuffs	160x100
40001-0000	1 PC	Abdominal Drape	150x150
CR 143038-003	1 PC	Maternity Pad	
CR 143014-000	1 PC	Cord Clamp Lock	
CR 141127-000	1 PC	Vaginal Packing Sponge	
CR 143121-000	1 PC	Bulb Syringe	
CR 143024-002	10 PCS	Gauze Swab X-Ray Detectable	10x10
CP 141081-004	1 PC	Placenta Basin	
CP 141066-000	1 PC	Carry Bag	

Delivery Pack II : Product Code - #66551-0002

MATERNITY MAT WITH POUCH





Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
CR 143005-000	2 PC	Apron	XL
CC 142023-000	1 PC	Maternity Mat with Pouch	
		• Under Pad	
		• Laminated Fabric	
CC 142005-000	2 PCS	Leggings	160x100
065-100-075S	1 PC	Absorbent Sheet	100x75
CR 143041-012	2 PCS	Mopping Pad	30x30
CR 143029-000	4 PCS	Hand Towels	30x30
065-075-075S	1 PC	Absorbent Baby Sheet	75x75
CR 143033-003	1 PC	ID Band Infant Pink Plain	
CR 143014-000	1 PC	Cord Clamp Lock	
CP 141066-000	1 PC	Carry Bag	



Delivery Pack III : Product Code - #66551-0003


MATERNITY MAT WITH FLUID COLLECTION POUCH




Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
CR 143027-024	2 PCS	Surgical Gowns	L
CR 143029-000	2 PCS	Hand Towels	
CR 143038-003	1 PC	Maternity Pad	
CC 142005-000	2 PCS	Leggings	160x100
CR 143038-004	1 PC	Maternity Mat with Fluid Collection Pouch	120x75
		• Maternity Absorbent Pad	
		• Fluid Collection Pouch	
		• Adhesive	
065-075-075S	2 PCS	Absorbent Sheet	75x75
CC 142013-000	1 PC	Absorbent Mat	30x30
CR 143005-000	2 PCS	Apron	XL
CR 143033-003	1 PC	ID Band Infant Pink Plain	
CR 143014-000	1 PC	Cord Clamp Lock	
CR 143084-000	1 PC	Under Pad	60x60
CR 143041-012	4 PCS	Mopping Pad	30x30
CP 141066-000	1 PC	Carry Bag	

Lithotomy Pack I : Product Code - #66563-0001

LITHOTOMY T^m SHEET



UNDER BUTTOCKS DRAPE WITH ABSORBENT PAD




Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
CR 143027-024	2 PCS	Surgical Gowns	L
CR 143029-000	2 PCS	Hand Towels	30x30
043-060-040A	2 PCS	Utility Drape	60x40
		• Adhesive	
CC 142005-000	2 PCS	Leggings	160x100
CR 143029-000	1 PC	Hand Towel	30x30
30226-0000	1PC	Lithotomy T ^m Sheet	190x150
		• 7x15 Cm Perineal Fenestration	
		• Adhesive	
		• Absorbent Laminated Reinforcement	
		• Tube Holder	
CR 143038-006	1 PC	Under Buttocks Drape with Absorbent Pad	
		• Absorbent Pad	
CP 141066-000	1 PC	Carry Bag	




Lithotomy Pack II : Product Code - #66563-0002


UNDER BUTTOCKS DRAPE WITH ABSORBENT PAD




SURGICAL SHEET



LEGGINGS



UTILITY DRAPE



Item Code	Qty.	Contents	Size(cm)
085-160-100S	1 PC	Back Table Cover Reinforced	160x100
CR 143027-024	1 PC	Surgical Gowns	L
CR 143029 - 000	2 PCS	Hand Towels	30x30
043-060-040A	2 PCS	Utility Drape	60x40
CR 143038-005	1 PC	Under Buttocks Drape with Absorbent Pad	
		• Absorbent Pad	
CC 142005-000	2 PCS	Leggings	160x100
043-160-1005	1 PC	Surgical Sheet	160x100
CP 141066-000	1 PC	Carry Bag	

